Quantachrome® ASiQwin™- Automated Gas Sorption Data **Acquisition and Reduction** © 1994-2017, Quantachrome Instruments version 5.21



9mm

<u>Analysis</u> Report

Operator: **Anton Paar** Date:2022/01/07 Operator: Chris Date:2022/01/24

Sample ID: Filename: CSH_GP01A_3_st1_2022_01_07.qps Comment: Sample Desc:

Sample Weight: 0.0337 g Autosorb iQ Station 1 Instrument:

Outgas Time: 0.0 hrs Outgas Temp.: 0°C Analysis gas: 6.58e-05 1/Torr Nitrogen Non-ideality: CellType:

Analysis Time: 14:24 hr:min Bath temp.: 77.35 K

Analysis Mode: Standard VoidVol Remeasure:off VoidVol. Mode: He Measure Cold Zone V: 5.20518 cc Warm Zone V: 7.69438 cc

dV(d)

Data Reduction Parameters

Thermal Transpiration: on Eff. mol. diameter (D): 3.54 Å Eff. cell stem diam. (d): 4.0000 mm


Calc. Model: N2 at 77 K on carbon (cylindr./sphere pores, QSDFT adsorption branch)

DFT method Rel. press. range: 0.0000 - 1.0000 Moving pt. avg: off

Adsorbate model Nitrogen Temperature 77.350K

Molec. Wt.: 28.013 **Cross Section:** 16.200 Å² **Liquid Density:** 0.806 g/cc





Pore width (Å)

DFT method summary

Pore volume = 0.154 cc/g

268.818 m²/g Surface area =

Lower confidence limit = 8.973 Å Fitting error =

0.969 %

Pore width (Mode) = 11.031 Å

Moving point average: off